

Tech Tip 0001 - *Reference Standards Product Grades*



Choosing the correct reference standard is essential to the success of your quality management program or research project. ChromaDex offers three different product grades to suit your technical and budgetary needs. These product grades are differentiated by the level of characterization and testing performed and reported on the certificate of analysis. The below information will guide you in choosing the right product grade for your use. If you have any questions about our product grades or would like assistance in your choice, please let us know.

Primary Analytical Standards (P)

Primary Standards are orthogonally characterized by multiple analytical techniques. These products include a complete certificate of analysis for each lot that lists adjusted (as-is) purity calculated by HPLC or GC chromatographic purity, Karl Fischer (water content), and GC (residual solvent). The identity of these standards is confirmed with Mass Spectrometry and NMR. The certificate of analysis includes analytical conditions, chromatograms, and spectra for your reference, allowing you to replicate the tests in your laboratory. When utilizing a reference standard for quantitative analysis, the best choice is a Primary Standard.

Analytical Standards (AS)

Analytical Standard grade is a step below our Primary grade regarding the level of characterization reported on the certificate of analysis accompanying each lot. The purity of these products is determined by HPLC or GC chromatographic purity. Identity is confirmed by one technique, typically Mass Spectrometry. Water and residual solvent content are not reported for Analytical Standards.

Reagent Grade Chemicals (RG)

Reagent Grade Chemicals are suitable as fine chemicals for many laboratories. However, they are not intended for use as reference standards. The certificate of analysis for RG chemicals contains basic physical properties, but the materials are not characterized for purity. As a result, RG chemicals are not suitable for direct quantitative use.

Tech Tip 0001 - Reference Standards Product Grades

ChromaDex product grades are differentiated by the level of characterization and testing performed and reported for each lot of material.

Grade	Chemical Data	Adjusted Purity	HPLC/GC %	Water %	Solvent %	NMR	Mass Spec
Primary Analytical Standard (P)	✓	✓	✓	✓	✓	✓	✓
Analytical Standard* (AS)	✓		✓				✓
Reagent Grade (RG)	✓						

Depending on your project and needs, different end uses require different product grades. Here is an overview of some typical uses and appropriate ChromaDex product grades for each.

Grade	Quantitative HPLC/GC	Quality Control	Method Devel.	Working Stnd.	Cell Assays	Basic Research
Primary Analytical Standard (P)	✓	✓	✓	✓	✓	✓
Analytical Standard (AS)	✓*	✓*	✓	✓	✓	✓
Reagent Grade (RG)					✓	✓

*Since Analytical Standard (AS) grade products are not corrected for water and solvent content, quantitative results from these standards will give content on chromatographic purity only and will not account for the water and solvent that is contained in the sample being tested. Whenever possible for the most accurate analytical results, the best choice for quantitative analysis is a Primary Analytical Standard (P).

In addition to our three product grades, ChromaDex offers two additional specialized types of Primary Standards. These products are full Primary Standards with additional value-added testing and results.

Primary Certified Reference Material (CRM)


For increased compliance if your laboratory is ISO accredited, ChromaDex produces select Primary Standards in accordance with our ISO 17034 Certified Reference Material Producer and ISO/IEC 17025 Chemical Testing Laboratory accreditations. Our Primary CRMs are the most rigorously tested materials in our catalogue, providing the best accuracy and most thorough certificates of analysis. In addition to the full Primary characterization, these materials include a reported uncertainty value for the HPLC or GC chromatographic purity. This overall value incorporates the uncertainty of the test method analysis, sample homogeneity, stability, and transport. The extra reported information will allow your laboratory to properly calculate measurement uncertainty for your sample analyses.

Compendial Traceable Primary Analytical Standards (CT)

As an added value to our clients, ChromaDex is now offering Primary Analytical Standards (P) that are also traceable to a compendial reference material. These products undergo full Primary characterization and also have a purity assay value tested against a compendial material such as supplied by the United States Pharmacopeia (USP), National Institute of Standards and Technology (NIST), or other similar international compendial. Certificates of Analysis for these products include full data for Primary Analytical Standards along with a listing of the compendial source, lot number, and assay value.

Tech Tip 0001 - Reference Standards Product Grades

Example: Primary Grade Certificate of Analysis



CONFIDENTIAL

Certificate of Analysis

For Laboratory Use Only

PRODUCT NAME Aloin B

PART NUMBER 00001626

STANDARD TYPE Primary (P)

LOT NUMBER 00001626-573

REPORT NUMBER CDXA-RSS-9170-00

SAMPLE NUMBER CDXP-20-00733

DATE OF SAMPLE 09/23/2020

DATE OF REPORT 10/21/2020

CHEMICAL NAMES (10R)-10- β -D-glucopyranosyl-1,8-dihydroxy-3-(hydroxymethyl)-9(10H)-anthracenone; Isobartalolin

CHEMICAL FORMULA C₂₁H₂₀O₈


MOLECULAR WEIGHT (MW) 418.39

CHEMICAL FAMILY Quinones

CAS NUMBER [28377-16-6]

EC(INECS) N/A

RECS CB1543000



ANALYTICAL RESULTS

TEST	METHOD	SPECIFICATION	RESULT
Adjusted Purity	NA	$\geq 85\%$	85.9%
HPLC	0.700.10.2.MERH22	$\geq 90\%$	96.6%
NMR	0.700.12.4	Conforms to structure	Conforms to structure
Mass Spectrum	0.700.12.2	Conforms	Conforms
Residual Solvent*	USP <467>	NA	Acetonitrile - 827 ppm 1,2-Dichloroethane - 93,500 ppm
Water	0.700.12.37	NA	1.8%
Appearance	NA	White powder	

*Testing performed at an ISO 17025 accredited subcontracted laboratory.

ADJUSTED PURITY: 85.9% IS BASED ON (100 - 9.3 SOLVENTS - 1.8 WATER) X 96.6% HPLC

STORAGE CONDITIONS


STORAGE -20 °C in a dry place.

EXPIRATION DATE 09/2025 under the above conditions.

This document is the property of ChromaDex, Inc. and/or its relevant affiliates and contains confidential and proprietary material for the sole use of the intended recipient(s). Any review, use, distribution or disclosure by or to others is strictly prohibited.

10900 Wilshire Blvd., Suite 600, Los Angeles, CA 90024 | T: +1 949-419-0288 | F: +1 949-419-0294 | www.chromadex.com
Laboratory address: 1751 S. Fordham St., Suite 300, Longmont, CO 80503

Page 1 of 6



CONFIDENTIAL

Certificate of Analysis

For Laboratory Use Only

ANALYTICAL CONDITIONS

INSTRUMENT AGILENT 1260 HPLC UV-VIS (DAD); AB-SCIEX 6400 TRIPLE TOF

COLUMN Phenomenex Luna C18(2) 150 x 2.00 mm, 5 μ m particle size

MOBILE PHASE A - 0.1% trifluoroacetic acid in Ultrapure water, B - Acetonitrile; Isocratic: 5% B for 5 minutes, then increasing to 95% B over 20 minutes. Hold at 95% B for 15 minutes, then decreasing to 5% B over 1 minute.

COLUMN TEMP. 60 °C


FLOW RATE 0.4 mL/minute

INJECTION VOL. 0.4 μ L

INJECTION CONC. 1.3 mg/mL in Methanol

DETECTION 205 \pm 4 nm


HPLC CHROMATOGRAM OF ALOIN B (CDXP-20-00733)



This document is the property of ChromaDex, Inc. and/or its relevant affiliates and contains confidential and proprietary material for the sole use of the intended recipient(s). Any review, use, distribution or disclosure by or to others is strictly prohibited.

10900 Wilshire Blvd., Suite 600, Los Angeles, CA 90024 | T: +1 949-419-0288 | F: +1 949-419-0294 | www.chromadex.com
Laboratory address: 1751 S. Fordham St., Suite 300, Longmont, CO 80503

Page 2 of 6




CONFIDENTIAL


Certificate of Analysis

For Laboratory Use Only

UV SPECTRUM AT RT = 13.1 MINUTES




HPLC CHROMATOGRAM OF DILUENT BLANK



This document is the property of ChromaDex, Inc. and/or its relevant affiliates and contains confidential and proprietary material for the sole use of the intended recipient(s). Any review, use, distribution or disclosure by or to others is strictly prohibited.

10900 Wilshire Blvd., Suite 600, Los Angeles, CA 90024 | T: +1 949-419-0288 | F: +1 949-419-0294 | www.chromadex.com
Laboratory address: 1751 S. Fordham St., Suite 300, Longmont, CO 80503

Page 3 of 6




CONFIDENTIAL


Certificate of Analysis

For Laboratory Use Only

MASS SPECTRUM OF ALOIN B (CDXP-20-00733)




¹H-NMR SPECTRUM OF ALOIN B (CDXP-20-00733) IN CD₃OD



This document is the property of ChromaDex, Inc. and/or its relevant affiliates and contains confidential and proprietary material for the sole use of the intended recipient(s). Any review, use, distribution or disclosure by or to others is strictly prohibited.

10900 Wilshire Blvd., Suite 600, Los Angeles, CA 90024 | T: +1 949-419-0288 | F: +1 949-419-0294 | www.chromadex.com
Laboratory address: 1751 S. Fordham St., Suite 300, Longmont, CO 80503

Page 4 of 6




CONFIDENTIAL

Certificate of Analysis

For Laboratory Use Only


¹³C-NMR SPECTRUM ALOIN B (CDXP-20-00733) IN CD₃OD



This document is the property of ChromaDex, Inc. and/or its relevant affiliates and contains confidential and proprietary material for the sole use of the intended recipient(s). Any review, use, distribution or disclosure by or to others is strictly prohibited.

10900 Wilshire Blvd., Suite 600, Los Angeles, CA 90024 | T: +1 949-419-0288 | F: +1 949-419-0294 | www.chromadex.com
Laboratory address: 1751 S. Fordham St., Suite 300, Longmont, CO 80503

Page 5 of 6




CONFIDENTIAL

Certificate of Analysis

For Laboratory Use Only

COSY-NMR SPECTRUM OF ALOIN B (CDXP-20-00733) IN CD₃OD



REVISION HISTORY

Revision History	Date of Revision	Document Changes
00	10/21/2020	New report


This document is the property of ChromaDex, Inc. and/or its relevant affiliates and contains confidential and proprietary material for the sole use of the intended recipient(s). Any review, use, distribution or disclosure by or to others is strictly prohibited.

10900 Wilshire Blvd., Suite 600, Los Angeles, CA 90024 | T: +1 949-419-0288 | F: +1 949-419-0294 | www.chromadex.com
Laboratory address: 1751 S. Fordham St., Suite 300, Longmont, CO 80503

Page 6 of 6

Tech Tip 0001 - Reference Standards Product Grades

Example: Analytical Standard Grade Certificate of Analysis




CONFIDENTIAL

Certificate of AnalysisFor Laboratory Use Only

PRODUCT NAME Valeric acid
PART NUMBER 00022195
STANDARD TYPE Analytical Standard (AS)
LOT NUMBER 00022195-WKF
REPORT NUMBER CDXA-RSS-7454-01
SAMPLE NUMBER CDXA-17-000507
DATE OF SAMPLE 01/11/2017
DATE OF RE-EVALUATION 03/05/2021
DATE OF REPORT 03/31/2021


CHEMICAL NAMES 1-Butanecarboxylic acid; n-Pentanoic acid; Pentanoic acid; Propylacetic acid
CHEMICAL FORMULA C₅H₁₀O₂
MOLECULAR WEIGHT (MW) 102.13
CHEMICAL FAMILY Oils
CAS NUMBER [109-52-4]
EC(ENECs) 203-677-2
ITECS YV610000

Structure



TEST	METHOD	SPECIFICATION	RESULT
GC	0.700.10.2.MERH1.40	2.75%	99.2%
Mass Spectrum	0.700.10.2.MERH1.40	Confirms	Confirms
Appearance	NA	NA	Liquid

STORAGE CONDITIONS
STORAGE +4 °C in a dry place.
EXPIRATION DATE 03/2026 under the above conditions.




Digitally signed by Elizabeth Hughes
Date: 2021.03.31 16:01:26 -06'00'

This document is the property of ChromaDex, Inc. and/or its relevant affiliates and contains confidential and proprietary material for the sole use of the intended recipient(s). Any review, use, distribution or disclosure by or to others is strictly prohibited.

10900 Wilshire Blvd., Suite 600, Los Angeles, CA 90024 | T: +1 949-419-0288 | F: +1 949-419-0294 | www.chromadex.com
Laboratory address: 1751 S. Fordham St., Suite 300, Longmont, CO 80503

Page 1 of 4




CONFIDENTIAL

Certificate of AnalysisFor Laboratory Use Only

ANALYTICAL CONDITIONS
INSTRUMENT AGILENT 7890B WITH FLAME IONIZATION DETECTOR (FID) AND AGILENT 8977B MASS SELECTIVE DETECTOR (MSD)
COLUMN Phenomenex Zebron ZB-Splus, 30 m x 0.32 mm x 0.25 µm
CARRIER GAS Helium
COLUMN TEMP. Isothermic 40 °C for 5.0 minutes, then increasing at 10.0 °C/minute to 340 °C. Hold at 340 °C for 5.0 minutes.
FLOW RATE 2.0 mL/minute
INJECTOR TEMP. 300 °C
INJECTOR SPLIT RATIO 20:1
INJECTION VOL. 1.0 µL
INJECTION CONC. 1.9 mg/mL in Chloroform
DETECTION FID (350 °C) & MSD (Transfer line 280 °C, Source 230 °C, Quad 150 °C)

GC CHROMATOGRAM OF VALERIC ACID (CDXA-17-000507)




This document is the property of ChromaDex, Inc. and/or its relevant affiliates and contains confidential and proprietary material for the sole use of the intended recipient(s). Any review, use, distribution or disclosure by or to others is strictly prohibited.

10900 Wilshire Blvd., Suite 600, Los Angeles, CA 90024 | T: +1 949-419-0288 | F: +1 949-419-0294 | www.chromadex.com
Laboratory address: 1751 S. Fordham St., Suite 300, Longmont, CO 80503

Page 2 of 4

Example: Reagent Grade Certificate of Analysis




CONFIDENTIAL

Certificate of AnalysisFor Laboratory Use Only


PRODUCT NAME Hederagenin
PART NUMBER 00008067
STANDARD TYPE Reagent grade (RG)
LOT NUMBER 00008067-815
REPORT NUMBER CDXA-RSS-7202-01
SAMPLE NUMBER CDXP-21-00204
DATE OF SAMPLE 07/13/2016
DATE OF RE-EVALUATION 03/25/2021
DATE OF REPORT 04/01/2021

CHEMICAL NAMES [36,4a]-3,23-Dihydroxyolean-12-en-28-oic acid; Astragaligenin E; Caulisapogenin; Hederagenic acid; Hederagenin
CHEMICAL FORMULA C₃₀H₄₈O₅
MOLECULAR WEIGHT (MW) 472.70
CHEMICAL FAMILY Triterpenes
CAS NUMBER [465-99-6]
EC(ENECs) 207-369-9
ITECS RK0177850

Structure



STORAGE CONDITIONS
STORAGE Room Temperature in a dry place.
EXPIRATION DATE 03/2026 under the above conditions.



Digitally signed by Elizabeth Hughes
Date: 2021.04.01 16:15:57 -06'00'

REVISION HISTORY


Revision History	Date of Revision	Document/Changes
00	08/03/2016	New report
01	04/01/2021	Passed re-evaluation by MS. Updated expiration date and removed melting point.

Note – Reagent Grade (RG) chemicals are not guaranteed as quantitative standards. This product line has been developed for research and qualitative purposes only.

This document is the property of ChromaDex, Inc. and/or its relevant affiliates and contains confidential and proprietary material for the sole use of the intended recipient(s). Any review, use, distribution or disclosure by or to others is strictly prohibited.

10900 Wilshire Blvd., Suite 600, Los Angeles, CA 90024 | T: +1 949-419-0288 | F: +1 949-419-0294 | www.chromadex.com
Laboratory address: 1751 S. Fordham St., Suite 300, Longmont, CO 80503


Page 1 of 1




CONFIDENTIAL

Certificate of AnalysisFor Laboratory Use Only

GC CHROMATOGRAM OF DILUENT BLANK




MASS SPECTRUM OF VALERIC ACID (CDXA-17-000507)



This document is the property of ChromaDex, Inc. and/or its relevant affiliates and contains confidential and proprietary material for the sole use of the intended recipient(s). Any review, use, distribution or disclosure by or to others is strictly prohibited.

10900 Wilshire Blvd., Suite 600, Los Angeles, CA 90024 | T: +1 949-419-0288 | F: +1 949-419-0294 | www.chromadex.com
Laboratory address: 1751 S. Fordham St., Suite 300, Longmont, CO 80503

Page 3 of 4



CONFIDENTIAL

Certificate of AnalysisFor Laboratory Use Only

REVISION HISTORY

Revision History	Date of Revision	Document/Changes
00	03/09/2017	New report
01	03/31/2021	Material passed reveal by GC-MS. Updated analytical results, analytical conditions, chromatograms, spectrum, and expiration date. Removed melting point.

This document is the property of ChromaDex, Inc. and/or its relevant affiliates and contains confidential and proprietary material for the sole use of the intended recipient(s). Any review, use, distribution or disclosure by or to others is strictly prohibited.

10900 Wilshire Blvd., Suite 600, Los Angeles, CA 90024 | T: +1 949-419-0288 | F: +1 949-419-0294 | www.chromadex.com
Laboratory address: 1751 S. Fordham St., Suite 300, Longmont, CO 80503

Page 4 of 4