



# **Best-in-Class Resins for Process Chromatography of Antibody**

**Tosoh Corporation**

**December, 2019**






# Content

- Market information in 2010 - 2024
- Downstream purification process for mAb
- Best-in-Class resins
- TOYOPEARL® properties
- Protein A, Protein L affinity chromatography
- Salt tolerant ion-exchange chromatography
- Hydroxyapatite (mixed-mode) chromatography
- Viral clearance
- Other mAb application
- Summary



# Worldwide Pharma Sales in Biotech vs. Traditional Technology in 2018-2024

Source; EvaluatePharma® May 2018



- Annual growth rate (CAGR) from 2018 to 2024 is 6.4 %.
- Pharma sales by biotech reaches 31 % in 2024.
- Pharma sales by biotech reaches 52 % in 2024 in top 100 products.
- Expected 12 mAbs become hit sales in top 20 pharma in 2024.



TOSOH

# Examples of New mAb Purification Process

Traditional

No salt HIC

Mixed Mode

Two step

Replace Pro A

Multiple Flowthrough

Genentech-2







Biogen (2013)

KBI Pharma (2016)

BTI (2016)

KIT (2018)

Kyowa Kirin (2017)

Astellas(2018)  
Millipore/Sigma

\* High load HIC

TOYOPEARL®

Hexyl-650C

\* Capto™ MMC

Nuvia™ cPrime

TOYOPEARAL

MX-Trp-650M

Sulfate-650F

(salt-tolerant)

\* MabSelect™ SuRe LX

Capto Adhere

\* TOYOPEARL

AF-rProtein A HC-650F

NH2-750F (salt-tolerant)

\* TOYOPEARL

GigaCap® S-650M

\* Depth Filter

Millistak® +A1HC,

Zeta Plus™ EXT, AC

\* TOYOPEARL

NH2-750F (salt-tolerant)

(\*Please remind any application may be patented before being public information.)



# Properties of TOYOPEARL Resins

- **Polymethacrylate polymer resin**
  - Pore size control ; MW range; 500 – 20,000,000
  - Particle size control; 30 – 150 um (grade S, F, M, C, EC)
  - High flow rate, large column packing
- **High capacity**
  - High capacity resin has an impact to decrease COGS
  - Grafting technology; IEC GigaCap series and other resins
  - IEC, HIC, MXC, AFC
- **Smaller particle grade**
  - TSKgel PW bulk IEC, HIC (20 um, 30 um)
  - Polishing step for protein, peptide and oligonucleotide
- **Prepacked columns for process development**
  - ToyoScreen® column (1, 5 mL), MiniChrom column (5 mL) for LC
  - RoboColumn® (200 uL, 600 uL) for Tecan microplate format
- **Supply chain**
  - Secure for large amounts supply by the 3rd manufacturing plants in 2019



# The Best-In-Class Resins on TOYOPEARL (High capacity, selectivity, purity and stability)

## Affinity chromatography

- **TOYOPEARL AF-rProtein A HC-650F**
  - High capacity, high alkaline stability, applicable to continuous chromatography
  - Large process column packing for multiple customers
- **TOYOPEARL AF-rProtein L-650F**
  - High capacity, alkaline durability for antibody fragment and scFv
  - High selectivity for bi-specific antibody

## Mixed-mode or salt-tolerant ion-exchange chromatography

- **TOYOPEARL NH<sub>2</sub>-750F**
  - High capacity, salt tolerant anion-exchanger (as mixed-mode)
  - High efficiency for removal of aggregate, endotoxin and viral clearance
- **TOYOPEARL Sulfate-650F**
  - High capacity salt tolerant cation-exchanger (as mixed-mode) with different selectivity
  - High efficiency for removal of aggregate, HCP and host DNA
- **Hydroxyapatite Ca<sup>++</sup>Pure-HA<sup>®</sup>**
  - Mixed-mode (Multi-mode) chromatography
  - Cost effective and easy transition from conventional hydroxyapatite



TOSOH

# Best-In-Class Resin Specification

Mode	TOYOPEARL			
	Affinity		Anion-exchanger	Cation-exchanger
Product	AF-rProtein A HC-650F	AF-rProtein L-650F	NH <sub>2</sub> -750F	Sulfate-650F
Base resin	Hydrophilic polymer			
Particle size	45 µm	45 µm	45 µm	45 µm
Ligand	Recombinant Protein A	Recombinant Protein L	Poly-amine	Sulfate
Ligand density	N/A	N/A	0.07~0.13 eq/L	≥ 0.53 eq/L
Static Binding Capacity (SBC; human IgG)	≥ 68 g/L-resin	≥ 64 g/L-resin	≥ 70 g/L-resin	≥ 114 g/L-resin
Storage	20 % Ethanol*			20 % Ethanol 0.2 mol/L sodium acetate
Storage temperature	2~8 degree C		Room temperature	

N/A; not applicable

\*\* For TOYOPEARL NH<sub>2</sub>-750F, washing with 0.1 mol/L HCl to substitute counter ions on the resin before storage.

	HYDROXYAPATITE	TOYOPEARL
	Mixed-mode (Multi-mode)	
Mode	Ca <sup>++</sup> Pure-HA	MX-Trp-650M
Base resin	Hydroxyapatite	Hydrophilic polymer
Particle size	39 µm	65 µm
Ligand	Ca <sup>2+</sup> , PO <sub>4</sub> <sup>3-</sup> , OH	Tryptophan
Ligand density	N/A	N/A
Static Binding Capacity (SBC; human IgG)	≥ 20 g/L-resin*	≥ 75 g/L-resin
Storage	Dry (or 1 mol/L NaOH)	20 % Ethanol
Storage temperature	Room temperature	

\* 55 g/L at 5 min residence time



# **Protein A & Protein L Affinity Chromatography**




# Binding of Protein A/Protein L to Antibody

## Protein A; Binding to Fc region

- Intact antibody (Tosoh protein A resin can also weakly bind to Fab region)
- Fc-fusion protein

## Protein L; Binding to kappa type light chain

- Antigen-binding fragment (Fab)
- Single-chain variable fragment (scFv)
- Single-domain antibody fragment (dAb)





# **TOYOPEARL AF-rProtein A HC-650F**



# TOYOPEARL AF-rProtein A HC-650F

- High capacity type
  - The highest binding capacity (similar to that on GE Prisma)
  - Higher (double) binding capacity at low temperature in a cold room against commercial Protein A resin
  - Applicable to weak binding antibody (e.g. IgM, Fab, Fc-fusion protein, mouse IgG)
- Purity
  - Improvement of removal of impurities with intermediate wash
  - Intermediate wash buffer at higher pH, addition of arginine or alcohol
- Alkaline stability
  - CIP over 300 cycles with 0.2 mol/L NaOH
- Installment of large process columns
  - 60 cm I.D. column ( over 50 L resin), secure supply chain
  - 15 cm height column gives better productivity



# Protein A Resin;

## Effective Intermediate Wash for Purification for mAb IgG<sub>2</sub>

Intermediate wash conditions	Sample Loading (mg/mL)	Yield (%)	HCP (ppm)	HMW (%)	Productivity (%)
Target (similar level to MabSelect SuRe LX)	35 - 49	> 85	< 1,400	< 3.0	-
Original; 25 mM Tris-HCl (pH 7.4) + NaCl (conductivity 12-14 mS/cm)	54	94	3,600	2.3	100
25 mM Tris-HCl (pH 9.0) + 10 % IPA	54	85	1,199	2.4	119
25 mM Tris-HCl (pH 9.0) + 0.5 mol/L arginine	54	89	951	2.5	125
25 mM Tris-HCl (pH 9.0)	54	93	1,725	2.5	131
25 mM Tris-HCl (pH 9.0) + 2 mol/L urea	54	81	1,118	2.5	114
25 mM Tris-HCl (pH 9.0) + 2 mol/L urea + 10 % IPA	54	65	439	2.2	91
25 mM Tris-HCl (pH 9.0) + 1 mol/L NaCl	54	84	1,953	2.5	119
25 mM Tris-HCl (pH 5.0) + 1 mol/L NaCl	54	81	2,415	2.4	113
25 mM Tris-HCl (pH 5.0) + 2 mol/L urea	54	74	1,429	2.4	104
25 mM Tris-HCl (pH 7.5) + 0.5 mol/L arginine	54	86	873	2.5	120

\* Ref.; K. Mehta et al., American Pharmaceutical Review, posted on Feb. 16, 2018

<https://www.americanpharmaceuticalreview.com/Featured-Articles/347357-Comparing-Performance-of-New-Protein-A-Resins-for-Monoclonal-Antibody-Purification/>

Data was modified from the original source.



# Dynamic Binding Capacity of mAb IgG<sub>1</sub> at Low Temperature on Protein A Affinity Chromatography

Temperature	Dynamic Binding Capacity (IgG <sub>1</sub> , g/L)			
	TOYOPEARL AF-rProtein A HC-650F		MabSelect SuRe	
	125 cm/hr	250 cm/hr	125 cm/hr	250 cm/hr
4 °C	33.2	27.8	18.6	11.0
22 °C	40.8	33.2	25.6	17.9
40 °C	40.2	34.2	30.8	24.7

Column size; 5 mm I.D. x 5 cm

Residence time; 2.4 min and 1.2 min, respectively

Data modified from the reference of W. Krepper et al., J. Chromatogr., A 1551 (2018) 59-68

- TOYOPEARL AF-rProtein A HC-650F showed the highest binding capacity at lower temperature in operation.



# TOYOPEARL AF-rProtein A HC-650F

## Chromatographic Conditions (1) Regular IgG

- **Binding and intermediate washing**

- 0.02 mol/L sodium phosphate buffer (pH 7.2), 0.15 mol/L NaCl
- When other impurities like green pigment is adsorbed, additional stepwise washing with buffer at pH above 4 may be effective.
- Intermediate washing with 50 mM Tris-HCl (pH 8.0), 1 mol/L NaCl \*<sup>1</sup>  
Or washing with 50 mM borate buffer (pH 9.0), 0.5 mol/L NaCl, 10 % isopropanol  
Or washing with 25 mM Tris-HCl (pH 7.5 - 9.0), 0.5 mol/L arginine, or 10 % isopropanol\*<sup>2</sup>
- Pre-washing with 0.1 mol/L NaCl (no buffer) would be also effective prior to low pH elution.

- **Elution**

- Elution pH is 0.2 point lower than other commercial Protein A resin
- 0.1 mol/L acetate or citrate (pH 3.7-3.0) or 0.1 mol/L glycine-HCl (pH 3.7-3.0)
- When aggregates were observed in eluted fraction, decrease concentration of acetate or citrate buffer below 0.06 mol/L.

- **Cleaning**

- Acidic solution of pH 2.7 - 2.5
- 0.1 mol/L NaOH, 0.5 mol/L NaOH is effective\*<sup>1</sup> to remove accumulated DNA/chromatin, while such cleaning conditions would deteriorate Protein A resin.
- For hydrophobic impurities, 20-30 % ethanol may be effective.

Reference \*<sup>1</sup>; P. Gagnon et al., J. Chromatogr., A, 1340 (2014)68-78

\*<sup>2</sup>; K. Mehta et al, American Pharmaceutical Review, posted in web on Feb. 16, 2018



# TOYOPEARL AF-rProtein A HC-650F

## Chromatographic Conditions (2) Fab, IgM, Mouse IgG

### ● Binding and intermediate washing

- Higher pH buffer without salt is recommended.
- 0.02 mol/L Tris-HCl (pH 8.4), no salt
- 0.02 mol/L sodium phosphate buffer (pH 7.2), no salt
- Pre-washing with 0.1 mol/L NaCl (no buffer) would be also effective prior to low pH elution.

### ● Elution

- Salt gradient from 0 to 0.5 mol/L NaCl in binding buffer
- 0.1 mol/L glycine-HCl (pH 3.7-3.0) or 0.1 mol/L acetate buffer (pH 3.7-3.0)
- When aggregates were observed in eluted fraction, decrease concentration of acetate or citrate buffer below 0.06 mol/L.

### ● Cleaning

- Acidic solution of pH 2.7-2.5
- 0.1 mol/L NaOH, to remove accumulated DNA/chromatin , 0.5 mol/L NaOH is effective<sup>\*1</sup>, while the cleaning conditions degradate Protein A resin.
- For hydrophobic impurities, 20-30 % ethanol may be effective.

<sup>\*1</sup>; Reference; P. Gagnon et al., J. Chromatogr., A, 1340 (2014)68-78

- Weak binding antibody and fragment would be applied to Protein L /Protein G affinity chromatography.



# **TOYOPEARL AF-rProtein L-650F**





# TOYOPEARL AF-rProtein L-650F

- High capacity type
  - The highest binding capacity
  - Higher (double) binding capacity at low temperature in a cold room against commercial Protein A resin
  - Artificial low molecule antibody (Fab, scFv, diabody)
- Purity
  - Improvement of removal of impurities by intermediate wash at higher pH, addition of arginine, alcohol, or denaturing agents
  - Intact bispecific antibody separation by both pH gradient and salt gradient at the same time
- Alkaline stability
  - Chemically stable with 0.1 mol/L NaOH for 10 hr
  - CIP over 300 cycles with 0.05 mol/L NaOH (combination with acidic pH 2.0 washing is effective.)
- Large scale purification
  - Process column for 10 L - 30 L at 15 cm height column



# Protein L Resin;

## Effective Intermediate Wash at pH 6.5 for Purification of Bispecific Tandem scFv from CHO Cell Culture Media



Loading amount (mg)	Wash pH	Additive for wash	Step elution pH	Recovery (%)	Monomer (relative %)	Aggregate (relative %)	Fragment (relative %)	HCP (ppm)	DNA (ppm)	Leached Protein L (ppm)*
Feed	-	-	-	-	24.6	63.2	12.3	14,782,068	2,676	-
5	6.5	-	2.85	86.9	94.5	5.1	0.4	7,341	7	<LOD
5	6.5	+ 0.25 M Arg	2.85	85.5	93.4	6.1	0.5	6,872	2	<LOD
5	6.5	+ 0.25 M Arg + 5% iPrOH	2.85	85.4	93.9	5.7	0.4	6,016	1	<LOD
20	6.5	-	2.85	83.4	83.4	16.1	0.4	10,522	3	<LOD
20	6.5	+ 0.25 M Arg	2.85	90.1	80.8	18.7	0.4	8,368	2	<LOD
20	6.5	+ 0.25 M Arg + 5% iPrOH	2.85	90.1	82.9	16.8	0.4	6,966	1	<LOD

Buffer; 0.1 mol/L sodium citrate buffer (pH 6.5)

LOD; Limit of detection


Column size; 1 mL

- By preliminary experiment, intermediate washing with 0.25 mol/L arginine was found effective to remove HCP from scFv fraction in CHO cell culture media.
- Additional isopropanol with arginine was found to be the most effective to remove HCP and to improve removal of other impurities.
- DBC of tandem scFv was calculated as 20 g/L by 10 % breakthrough curve at 0.5 g/L of the purified sample.
- In case of large sample loading, aggregate contents increased with decrease of monomer.




# SEC analysis of Tandem scFv purification

## Optimized Conditions (5 mg Load) in purification of Tandem scFv by TOYOPEARL AF-rProtein L-650F



Before purification



After AF-rProtein L-650F  
one step purification  
for 5 mg (20 mg) load

Tandem scFv monomer  
purity; 93.9 % (82.9%)



# TOYOPEARL AF-rProtein L-650F

## Chromatographic Conditions for Single Chain Fragment (scFv)

- **Adsorption and intermediate wash**

- (1) **scFv, tandem scFv (from CHO cell)**

- 0.02 - 0.1 mol/L sodium citrate buffer (pH 6.5)
    - 0.02 - 0.1 mol/L sodium citrate buffer (pH 6.5), 0.25 - 1.0 mol/L arginine + 5 - 10 % isopropanol

- (2) **scFv (inclusion body from *E. coli*)**

- 0.1 mol/L sodium phosphate buffer (pH 6.5)
      - 0.1 mol/L sodium phosphate buffer (pH 6.5), 2 mol/L guanidine/HCl
- (Reference; J. Vajda et al., Tosoh Bioscience Technical Note)

- **Elution**


- 0.1 mol/L citrate or glycine-HCl (pH 2.0 – pH 2.85)

- **Cleaning**

- Acidic solution described above at pH 2.0 – pH 2.2
  - CIP; 0.05 - 0.1 mol/L NaOH (limited alkaline washing; contact time for 15 min or 1 CV)
  - For cleaning of hydrophobic impurities; 10-20% ethanol or isopropanol



# CIP of TOYOPEARL AF-rProtein L-650F with 0.05 mol/L NaOH



- For not so alkaline stable resin, acidic wash like citric acid (pH 2.2-2.0) followed by 0.05 mol/L NaOH gives more effective cleaning than only 0.1 mol/L NaOH washing
- TOYOPEARL AF-rProtein L-650F is stable over 300 CIP cycles with 0.05 mol/L NaOH


Ref.; N. Eifler et al., Biotechnol. Prog., (2014) Vol. 30 No. 6, p.1311-1318



# **Salt-Tolerant Ion-Exchange Mixed-Mode (Multi-Mode) Chromatography**

- **TOYOPEARL NH<sub>2</sub>-750F**
- **TOYOPEARL Sulfate-650F**
- **Hydroxyapatite Ca<sup>++</sup>Pure-HA**
- **TOYOPEARL MX-Trp-650M**

# Global Market in Mixed (Multi)- Mode Chromatography Resin




Global market size; USD 75.4 M in 2017

USD 130.3 M by 2025 (CAGR;7%)

Ref.; Grand View Research; Mixed Mode Chromatography Resin Market Analysis Report By Technique (Ion exchange - Hydrophobic, Hydroxyapatite), By End Use (Pharmaceutical & Biotechnology), And Segment Forecasts, 2018 - 2025

<https://www.grandviewresearch.com/industry-analysis/mixed-mode-chromatography-resin-market>

# Salt-Tolerant Ion-Exchangers and Multi (Mixed) Mode Chromatography Resins

Resin	Functional group	Separation mode	Target molecule
TOYOPEARL NH <sub>2</sub> -750F	Polyamine -NH <sub>2</sub> (grafted)	<ul style="list-style-type: none"> <li>• Anion-exchange</li> <li>• Hydrophobic interaction</li> </ul>	Serum albumin
			Antibody (monomer, aggregate, scFv)
			Endotoxin, DNA (removal)
			Virus (removal)
TOYOPEARL Sulfate-650F	Sulfate -O-SO <sub>3</sub> (grafted)	<ul style="list-style-type: none"> <li>• Cation-exchange</li> <li>• Hydrophobic interaction</li> </ul>	Antibody (monomer, aggregate)
			small antibody (scFv)
			Antithrombin III
TOYOPEARL MX-Trp-650M	<b>Tryptophane</b> 	<ul style="list-style-type: none"> <li>• Cation-exchange</li> <li>• Hydrophobic interaction</li> </ul>	Antibody (IgG, Bispecific, aggregate)
			Small antibody (scFv, Fab, Fc)
			Serum protein, albumin (recombinant)
			Enzyme
			Insulin-like growth factor (IGF-1)
Ca <sup>++</sup> Pure-HA <sup>®</sup>	Hydroxyapatite Ca <sub>10</sub> (PO <sub>4</sub> ) <sub>6</sub> (OH) <sub>2</sub>	<ul style="list-style-type: none"> <li>• Cation-exchange</li> <li>• Metal chelate</li> </ul>	Antibody, aggregate
			Antibody fragment
			Enzyme





# Salt-Tolerant Anion-Exchanger (Mixed-Mode); TOYOPEARL NH<sub>2</sub>-750F

## Properties

- Salt-tolerant anion-exchanger with huge pore on resin, high capacity
- High selectivity with polyamine functional group (mixed-mode effect)
- Removal HMW impurities like aggregate, virus, endotoxin and DNA
- Flow-through purification after Protein A and/or cation-exchanger
- Application at various pH and salt concentration


## Hints

- Optimization of salt concentration and pH in buffer
  - Applicable to flow-through mode at over 0.15 mol/L NaCl in most cases
- Avoid multivalent anion buffer like phosphate, sulfate, carbonate, citrate, etc. for stronger binding
- Storage of resin with 20 % ethanol solution after washing with acidic solution at pH 4



# Optimizing Flow-Through (FT) Conditions on TOYOPEARL NH<sub>2</sub>-750F


- Load sample
  - Aggregate; 3.1%
  - CHO-HCP; 3,422 ppm
  - DNA; 0.70 %



Ref.; W. Evans, PITTCOON 2019, Philadelphia PA

TOSOH BIOSCIENCE

Conditions selected  
for testing



# Analysis of Flow-Through Fraction by SEC

	Recovery (monomer)	Ratio				
		HMW1 (Aggregate)	HMW2 (Dimer)	Monomer	LMW1 (Fragment)	LMW2 (Fragment)
Load		10.2 %	6.5 %	76.6 %	4.7 %	2.0 %
30 g/L	76.2 %	0.1 %	3.5 %	90.7 %	3.6 %	2.1 %
60 g/L	86.5 %	1.2 %	4.7 %	87.6 %	4.4 %	2.0 %
90 g/L	95.3 %	6.6 %	6.3 %	80.7 %	4.5 %	1.8 %
150 g/L	98.7 %	8.1 %	6.7 %	78.7 %	4.7 %	1.7 %

- Excellent removal of HMW1 impurity at  $\leq 45$  mg/mL load ratio, breakthrough at  $> 60$  mg/mL
- Significant HMW 1 breakthrough noted starting at 60 mg/mL load ratio
- Poor fragment removal

## Conditions

Column; TOYOPEARL NH<sub>2</sub>-750F, 5 mm I.D. x 5 cm (1 mL)

Flow through; 50 mmol/L MES-0.1 mol/L acetate (pH 6.5), at 9 mS/cm

Elution; 50 mmol/L MES + 1.0 mol/L NaCl (pH 6.5)

Sanitization: 0.5 mol/L NaOH

Flow rate; Flow rate: 0.25 mL/min (4 min residence time)

Temperature: ambient

Detection: UV (280 nm), conductivity (mS/cm)

Sample; mAb-01A, 7.5 mg/L, 9 mS/cm, pH 6.5, 20 mL (150 mg/mL load ratio)

Instrument: ÄKTA avant 25 (Unicorn 7.3)

Ref.; W. Evans, PITTCOON 2019, Philadelphia PA (data modified)



# Salt-Tolerant Cation-Exchanger (Mixed-Mode); TOYOPEARL Sulfate-650F

## Properties


- Salt-tolerant cation-exchanger with larger pore size, high capacity
- High selectivity by Sulfate functional group (mixed-mode effect)
- Best removal of aggregate antibody (HMW) in cation-exchangers
- Removal of impurities like virus, endotoxin and DNA
- Bind/Elute separation after Protein A and anion-exchanger
- Heparin-like affinity interaction with sulfate functional groups

## Hints

- Optimization of buffer pH and salt concentration
  - In most cases, buffer at pH 5.0 or lower, salt-tolerant separation can be achieved at 0.15 mol/L NaCl or higher.
- Operation at pH 4 - 12, CIP/short term cleaning at pH 3 - 13
- Storage of resin with 20 % ethanol including 0.2 mol/L sodium acetate at neutral or weak acidic pH



# Salt-Tolerant Cation-Exchanger (Mixed-Mode) TOYOPEARL Sulfate-650F



Column; TOYOPEARL Sulfate-650F  
(7.5 mm I.D. X 7.5 cm)

Elution;

Buffer A: 0.054 mol/L acetate buffer (pH 5.5)

Buffer B: 1.0 mol/L NaCl + Buffer A

Linear gradient:

Buffer A (95%)- Buffer B (100%), 58.5 min

Flow rate; 1.0 mL/min

Temperature; ambient


Sample; 3 g/L monoclonal humanized IgG  
acid/high temperature treatment  
(Injection volume 90  $\mu$ L)

- TOYOPEARL Sulfate-650F showed best resolution for antibody between monomer and aggregate.



# New Salt-Tolerant Cation-Exchanger; mAb Aggregate Removal by Bind/Elute Method

- Salt tolerant
- High capacity
- **Aggregate removal**
- Viral clearance
- Heparin-like affinity




Purification results	Protein A	Sulfate-650F
Sample loading (mg/mL resin)	48	101
Recovery, monomer (%)	n/a	83
Dimer (%)	3.90	2.40
Aggregate (%)	0.54	0.07
HCP (ppm)	1,260	134
Protein A (ppm)	1.20	0.04

\*Elution at 260 mmol/L NaCl

Resin: TOYOPEARL Sulfate-650F  
Column: 6.6 mm ID × 3.0 cm (1.0 mL)  
Mobile phase: A: 50 mmol/L acetate-Tris, 100 mmol/L NaCl, pH 5.2  
B: 50 mmol/L acetate-Tris, pH 5.2, NaCl as indicated  
C: 50 mmol/L acetate-Tris, 1.0 mol/L NaCl, pH 5.2  
D: 0.5 mol/L NaOH  
Flow rate: 45 cm/hr (4 min residence time)  
Detection: UV @ 280 nm (mAU)  
Temperature: ambient  
Injection vol.: 5.3 mL (97 mg/mL-resin load ratio)  
Sample: TBL-mAb-01, 19.1 mg/mL  
Instrument: ÄKTA avant 25

# Adsorption Capacity for scFv on Various Cation-Exchangers after Desalting



**Figure 4.** Screening of macropore ion exchange resins for equilibrium binding capacity of scFv. Hatched bars: direct capturing of scFv from refolding solution with anion exchangers; gray bars: capturing of scFv from desalted refolding solution with cation exchangers.

- **TOYOPEARL Sulfate-650F showed the best binding capacity for scFv.**

Ref. ; N. Walch et al., Biotechnol. J., 2017, 12, 1700082



# Multi (Mixed) -Mode (Cation-Exchange Type) Resin; TOYOPEARL MX-Trp-650M

## Properties

- Tryptophan (amino acid) immobilized
- Multi-mode (Cation-exchange/hydrophobic interaction) separation
- High capacity; IgG  $\geq$  75 g/L
- Sharper elution peak and higher recovery
- Salt-tolerant separation for antibody and other proteins


## Hint for separation

- Optimization of salt concentration and pH in eluent
- Adsorption at pH 4.0 - pH 5.5, elution at pH 7.0 - 8.0 with salt
- CIP, short term cleaning at pH 3 - pH 13
- Storage by shading container, ethanol solution without oxidative reagent





# Purification of Monoclonal Antibody from Cell Culture Media by TOYOPEARL MX-Trp-650M



## Conditions

Column; TOYOPEARL MX-Trp-650M (6 mm I.D. x 4 cm)

Binding buffer; 0.05 mol/L acetate buffer (pH 4.7) + 0.1 mol/L NaCl

Elution buffer; 0.1 mol/L Tris-HCl buffer (pH 8.5) + 0.3 mol/L NaCl

Flow rate (binding); 1.0 mL/min (212 cm/h)

(elution) ; 2.0 mL/min (at 124 min)


Detection; UV(280 nm)

Sample; CHO cell culture media, monoclonal antibody (1 mg/mL), diluted with binding buffer

Sample conductivity; 12 mS/cm

# Hydroxyapatite Ca<sup>++</sup>Pure-HA

- Media is composed of a unit cell comprised of  $\text{Ca}_{10}(\text{PO}_4)_6(\text{OH})_2$
- 10 nm x 100 nm hexagonal cross section crystals are agglomerated into particles and then heated to create stable welds at the crystal contact points
- Spherical, macroporous form of a hexagonal crystalline structure
- Sintered at high temperatures
- for increased mechanical and
- chemical stability






# Properties on Hydroxyapatite Ca<sup>++</sup>Pure-HA

- Effective to higher MW biomolecule like antibody by large 50 nm pore size of matrix
- **DBC on antibody; 55 g/L** (IgG; 2 g/L, residence time; 5 min)
- Similar or better capacity and mechanical stability than commercial hydroxyapatite
- CIP over 65 cycles with 1 mol/L NaOH  
(operational buffer pH above 6.5, unstable in acidic condition)
- **Mixed-mode matrix (cation-exchange type) effective to separate antibody and removal of impurities**
- **30 % lower price than commercial hydroxyapatite to replace existing process matrix**
- Target molecules  
Antibody and its fragment, recombinant protein, enzyme, virus, vaccine, nucleic acid (plasmid, double stranded DNA)

# Excellent Selectivity on Hydroxyapatite

## *Aggregate removal from purified IgG using **KCl** for elution*



- $\text{Ca}^{++}$ Pure-HA is compatible with potassium phosphate buffer as a loading buffer and potassium chloride as an elution buffer at different pH.
- Similar mAb aggregate removal performance across all three pH values.
- The separation profiles show high resolution between the monomer peak and the aggregate peaks, and increased pH makes earlier elution



# High Purity & Yield for mAb Purification

Salt	pH	Peak molarity (mmol/L)	Pool Volume (CV)	Recovery monomer (%)	Aggregate (%)	Fragment (%)	CHO-HCP (ppm)
Load					6.6	0.6	2,200
KCl	6.5	814	5.3	73	1.3	0.5	500
	7.0	615	3.9	80	1.8	0.3	590
	7.5	509	3.9	81	2.2	0.3	910

- > 80% recovery of mAb monomer at optimal pH (7.0)
- KCl elution buffer shows good aggregate removal across pH conditions.
- This elution buffer shows that the lower eluting buffer pH gives better aggregate removal.
- Ca<sup>++</sup>Pure-HA significantly reduces host cell proteins (HCP) at pH of 6.5 and 7.0.



# Attention fo Usage of Ca<sup>++</sup>Pure-HA

- Some limitation of mobile phase usage due to patent infringement
- Gradient elution by salt
  - Use 25 mM or higher concentration of phosphate buffer
  - Use KCl as salt (avoid use of NaCl)  
(KCl has better solubility than NaCl for salt gradient elution.)
- Gradient elution by phosphate buffer
  - No limitation for usage
- Neutralization of hydroxyapatite in initial cleaning
  - At first cleaning should be achieved with phosphate buffer, not to use NaOH solution for neutralization of hydronium ion in hydroxyapatite.



# Selection of TOYOPEARL Ion-Exchangers and/or Hydroxyapatite (HA)

## Customer already using hydroxyapatite

- Recommend **Ca<sup>++</sup>Pure-HA**
- **Cost advantage (30 % less price)**
- No significant difference in separation  
(replacement of existing matrix)
- Slightly better binding capacity
- Slightly better mechanical stability
- FOC sample demonstration  
(Anyway, ask trial resin with bulk matrix or MiniChrom)

## Customer not yet used hydroxyapatite

- Recommend **TOYOPEARL Sulfate-650F**, or **MX-Trp-650M**, (or **NH<sub>2</sub>-750F**, in case)
- TOYOPEARL is chemically stable in alkaline and acid.  
(HA is unstable in acidic solution.)  
(HA only is used at pH above 6.5)
- FOC sample demonstration  
(Anyway, ask trial resin with bulk matrix or MiniChrom)



# Summary of Viral Clearance of Antibody by Chromatographic Resins

Chromatography Mode	Resin	Separation	Log Reduction Value (LRV)*			
			MuLV	PRV	Reo	MVM
Anion-exchange (AIEC)	TOYOPEARL GigaCap Q-650M	Flow through	5.8	5.7	5.7	6.0
Cation-exchange (CIEC)	TOYOPEARL GigaCap S-650M	Bind-elute	4.3	5.5	1.4	-0.2
Mixed-mode (MXC)	TOYOPEARL NH <sub>2</sub> -750F**	Flow through	4.9	5.9	5.7	6.4
Mixed-mode (MXC)	TOYOPEARL Sulfate-650F**	Bind-elute	4.7	5.4	6.0	1.2
Mixed-mode (MX)	TOYOPEARL MX-Trp-650M	Bind-elute	2.8	3.5	1.8	0.6
Mixed-mode (MXC)	Hydroxyapatite Ca <sup>++</sup> Pure-HA	Flow through	3.3	2.3	2.7	-0.3
Hydrophobic interaction (HIC)	TOYOPEARL PPG-600M	Bind-elute	2.4	0.9	4.8	0.8
Hydrophobic interaction (HIC)	TOYOPEARL Hexyl-650C	Flow through	3.4	1.8	6.2	0.0
Affinity (AFC)	TOYOPEARL AF-rProtein A HC-650F	Bind-elute	1.7	1.2	0.0	2.2
Affinity (AFC)	TOYOPEARL AF-rProtein L-650F	Bind-elute	2.5	1.7	3.0	2.5

\* LRV was measured independently for each product by optimized conditions for model sample as monoclonal antibody or albumin (Hexyl-650C).

\*\*Salt-tolerant ion-exchanger

- Anion-exchanger and salt-tolerant (mixed-mode) resins show better viral clearance.






# Other New Application to mAb Separation

- Bispecific mAb separation
  - One-step purification by [TOYOPEARL AF-rProtein L-650F](#) by linear gradient of both pH and salt, (Chugai, 2018, PCT WO 2018/159615 A1)
  - Two-step purification by [TOYOPEARL MX-Trp-650M](#) and [Butyl-600M](#)  
(J. Vadjia et al., Bioprocess Int. 16 (10) (2018) 59-62), Novimmune SA, US patent 2016/0264685 A1)
- Antibody-Drug Conjugate (ADC) separation by HIC
  - Separation by less hydrophobic resin, [TOYOPEARL PPG-600M](#) or [Ether-650M](#)
- Continuous production of mAb
  - PEG treatment followed by [TOYOPEARL GigaCap S-650M](#)
  - [TOYOPEARL AF-rProtein A HC-650F](#) followed by AEC and MXC
- Flow-through HIC application
  - New flow-through HIC resin for mAb (trial resin)
  - Removal of detergent in cell culture feedstock prior to Protein A affinity purification on [TOYOPEARL Butyl-600M](#)




# New Process Development for Bispecific mAb

TOSOH







# Purification of Bispecific mAb by Multi (Mixed)-Mode and HIC



TOYOPEARL  
MX-Trp-650M




TOYOPEARL  
Butyl-600M




Ref.; J. Vadja et al., BioProcess International, October (2018) 16 (10) 59-62

# Purification of ADC by HIC on TOYOPEARL PPG-600M




Ref. V. Leur et al., Bioprocess International, 16(12)E2, Dec. 2018, modified



Column; ToyoScreen® PPG-600M 6.4 mm I.D. x 3 cm (1 mL)  
 Elution; A: 25 mmol/L sodium phosphate (pH 7.0)  
 + 1.5 mol/L ammonium sulfate  
 B; 25 mmol/L sodium phosphate (pH 7.0)/iso-propanol  
 = 75/25  
 60 min linear gradient from Buffer A to B  
 Flow rate; 1 mL/min, Detection; UV (280 nm)  
 Sample; ADC (1 g/L, 50  $\mu$ L)

# Separation of Antibody Drug Conjugate Mimic by HIC on TOYOPEARL PPG-600M



## Conditions

Column; TOYOPEARL PPG-600M (6.6 mm I.D. x 10 cm)

Elution; Buffer A: 0.1 mol/L sodium phosphate (pH 6.5) containing 1.5 mol/L ammonium sulfate

Buffer B: 0.1 mol/L sodium phosphate (pH 6.5)

Step gradient, 100 % Buffer A (5CV), 70 % Buffer B (5CV), 80 % Buffer B (5 CV), 100 % Buffer B (5 CV)

Flow rate; 250 cm/hr for equilibration and washing, 175 cm/hr for step elution

Detection; UV (280 nm), VIS (495 nm)

Sample; adalimumab, conjugated with FITC (5 mg/mL resin)

Ref.; Tosoh Bioscience Application Note (2019)



# Summary

- Best-in-Class process chromatography resins as of high capacity, high selectivity and alkaline stability for purification of antibody.
- TOYOPEARL Protein A and Protein L affinity chromatography, and salt-tolerant ion-exchangers showed good purification of various antibodies like intact IgG and single chain fragments (scFv).
- Combination of capture step, Protein A/Protein L affinity chromatography and salt-tolerant ion-exchangers would be powerful tool to purify monoclonal antibody and to remove impurities like aggregate, HCP, host DNA, virus and Protein A/Protein L leachable.



# Commercial Ion-Exchangers and Mixed-Mode Resins

Resin particle size	TOSOH	GE Life Science	Merck	Bio-Rad	Others
AEC 50-200 µm	TOYOPEARL GigaCap Q-650M TOYOPEARL QAE-550C TOYOPEARL Q-600C AR TOYOPEARL NH <sub>2</sub> -750F	Capto Q Q Sepharose XL Q Sepharose Fast Flow Q Sepharose Big Beads MacroCap Q CaptoCore Q	Eshumuno Q Fractogel EMDTMAE Hicap Fractogel EMD TMAE (M)	Nuvia Q MacroPrep High Q UNOsphere Q	POROS 50 micron HQ Pall Q Hypercel YMC BioPro QA
AEC 20-35 µm	TOYOPEARL GigaCap Q-650S TSKgel SuperQ-5PW (20), (30)	Capto Q ImPress Q Sepharose HP Source Q15, Q30	Fractogel EMD TMAE (S)	MacroPrep 25 Q	
AEC 50-100 µm	TOYOPEARL GigaCap DEAE-650M TOYOPEARL DEAE-650M/C	Capto DEAE DEAE Sepharose Fast Flow	Fractogel EMD DEAE	MacroPrep DEAE	
AEC 20-35 µm	TOYOPEARL DEAE-650S TOYOPEARL GigaCap DEAE-650S TSKgel DEAE-5PW (20), (30)	DEAE Sepharose HP			
CEC 50-200 µm	TOYOPEARL GigaCap S-650M TOYOPEARL SP-550C TOYOPEARL SP-650M/C TOYOPEARL MegaCap II SP-550EC TOYOPEARL Sulfate-650F	Capto S SP Sepharose XL SP Sepharose Fast Flow SP Sepharose Big Beads MacroCap SP	Eshumuno S Fractogel EMD SE Hicap (M) Fractogel EMD SO3-	Nuvia S MacroPrep High S UNOsphere S	POROS XS POROS 50 micron HS Pall S Hypercel YMC BioPro SP
CEC 20-35 µm	TOYOPEARL GigaCap S-650S TOYOPEARL SP-650S TSKgel SP-5PW (20), (30)	Capto S ImPress S Sepharose HP Source S15, S30	Fractogel EMD SO3-(S)	MacroPrep 25 S	
CEC 50-100 µm	TOYOPEARL GigaCap CM-650M TOYOPEARL CM-650M/C	CM Sepharose Fast Flow	Fractogel EMD COO-	MacroPrep CM	
CEC 20-35 µm	TOYOPEARL CM-650S	CM Sepharose HP			
MMC 50-100 µm	(TOYOPEARL NH <sub>2</sub> -750F)	Capto Adhere			Pall MEP HyperCel HEA HyperCel, PPA HyperCel
MMC 50-100 µm	TOYOPEARL MX-Trp-650M (TOYOPEARL Sulfate-650F)	CaptoMMC	Eshumuno HCX	Nuvia cPrime	
Hydroxyapatite 20-55 µm	Ca <sup>++</sup> Pure-HA			CHT Ceramic hydroxyapatite Type-I, II	



# BioSeparation on TSKgel® Columns and TOYOPEARL® Resins

Mode	Target	TSKgel Columns	TOYOPEARL, Preparative HPLC Resins	Remarks
Size-Exclusion Chromatography (SEC)	Peptide	TSKgel UP-SW2000 (2019) TSKgel SuperSW2000 TSKgel G2000SWxL	TOYOPEARL HW-40S/F	•High resolution on TSKgel SW series •Separation range on TOYOPEARL HW-40 similar to TSKgel G2500PWxL
	Protein	TSKgel UP-SW3000 TSKgel SuperSW mAb HR/HTP TSKgel G3000SWxL	TOYOPEARL HW-50S/F TOYOPEARL HW-55S/F	•High resolution on TSKgel SW series. •Separation range on TOYOPEARL HW-55 similar to TSKgel G3000PW
	Aggregate High Molecular Bio-complex	TSKgel UltraSW Aggregate, G4000SWxL TSKgel G4000PWxL, G5000PWxL TSKgel G6000PWxL/G-DNA-PW	TOYOPEARL HW-65S/F TOYOPEARL HW-75S/F	•TSKgel UltraSW Aggregate separates aggregates up to 5 mer •TOYOPEARL HW-75 has higher pore size.
	Oligosaccharide Polysaccharide	TSKgel SuperOligo-PW®, G-Oligo-PW TSKgel G2500PWxL, G3000PWxL TSKgel G5000PWxL, G6000PWxL TSKgel GMPWxL, SuperMultiporePW-H, -M	TOYOPEARL HW-40S/F TOYOPEARL HW-50S/F TOYOPEARL HW-55S/F TOYOPEARL HW-65S/F, HW-75S/F	•TSKgel Oligo-PW series for non-ionic saccharide •Separation range on TOYOPEARL HW-40 up to 8mer
Ion-Exchange Chromatography (IEC)	Acidic Molecule (Strong Anion-Exchange)	TSKgel Q-STAT, DNA-STAT, DNA-NPR® TSKgel SuperQ-5PW TSKgel BioAssist® Q	TOYOPEARL GigaCap® Q-650M/S TOYOPEARL SuperQ-650S/M TOYOPEARL Q-600C AR, QAE-550C TSKgel SuperQ-5PW(20), (30)	•TOYOPEARL GigaCap series with higher adsorption capacity •TOYOPEARL Q-600C AR with higher alkaline stability
	Acidic Molecule (Weak Anion-Exchange)	TSKgel DEAE-NPR TSKgel DEAE-5PW	TOYOPEARL GigaCap DEAE-650M TOYOPEARL DEAE-650S/M TOYOPEARL NH2-750F TSKgel DEAE-5PW (20), (30)	•TOYOPEARL GigaCap series with higher adsorption capacity •TOYOPEARL NH2-750F is salt-tolerant ion-exchanger (and mixed-mode) with higher resolution of huge bio-molecule
	Basic Molecule (Strong Cation-Exchange)	TSKgel SP-STAT TSKgel SP-NPR TSKgel SP-5PW	TOYOPEARL GigaCap S-650M TOYOPEARL SP-650S/M, SP-550C TOYOPEARL MegaCap II SP-550EC TOYOPEARL Sulfate-650F TSKgel SP-5PW(20), (30), SP-3PW(30)	•TOYOPEARL GigaCap series with higher adsorption capacity •TOYOPEARL Sulfate-650F is salt-tolerant ion-exchanger (mixed-mode). •TOYOPEARL SP-550 series and TSKgel SP-3PW(30) for insulin separation
	Basic Molecule (Weak Cation-Exchange)	TSKgel CM-STAT TSKgel CM-5PW	TOYOPEARL GigaCap CM-650M TOYOPEARL CM-650S/M	•TSKgel CM-STAT with higher resolution
Mixed-Mode Chromatography (MXC)	Acidic Molecule (Weak Anion-Exchange)	N/A	TOYOPEARL NH2-750F	•Salt-tolerant ion-exchanger
	Basic Molecule (Cation-Exchange) (Hydroxyapatite)	N/A	TOYOPEARL Sulfate-650F TOYOPEARL MX-Trp-650M Ca <sup>++</sup> Pure-HA	•TOYOPEARL Sulfate-650F is salt-tolerant ion-exchanger •Hydroxyapatite
Hydrophobic-Interaction Chromatography (HIC)	High Hydrophobic Protein	TSKgel Ether-5PW	TOYOPEARL Ether-650M, PPG-600M	•TOYOPEARL 600 series have higher binding capacity for antibody.
	Medium Hydrophobic Protein	TSKgel Butyl-NPR TSKgel Phenyl-5PW	TOYOPEARL Phenyl-600M, Phenyl-650S/M TOYOPEARL Butyl-600M, Butyl-650S/M TOYOPEARL SuperButyl-550C	•TOYOPEARL 600 series have higher binding capacity for antibody. •TOYOPEARL 550 series for low molecule protein separation
	Low Hydrophobic Protein Nucleic Acid	N/A	TOYOPEARL Hexyl-650C	•Plasmid DNA and oligosaccharide separation
Affinity Chromatography (AFC)	Antibody (Fc, glycan) Fc Fusion Protein	TSKgel FcR-IIIa-NPR TSKgel Protein A-5PW	TOYOPEARL AF-rProtein A HC-650F	•TSKgel FcR-IIIa-NPR can separate antibody according to activity.
	Antibody (Kappa Light)	N/A	TOYOPEARL AF-rProtein L-650F	•For artificial small antibody
	Heparin-Affinity Molecule	TSKgel Heparin-5PW	TOYOPEARL AF-Heparin HC-650M	•For plasma protein with heparin affinity
	Glycoprotein (Diol compound)	TSKgel Boronate-5PW	N/A	•For glycosylated protein separation
	Protein, Peptide	TSKgel Chelate-5PW	TOYOPEARL AF-Chelate-650M	•Metal ion chelated affinity chromatography (IMAC)
	Nucleotide-Dependent Affinity	N/A	TOYOPEARL AF-Red-650M	•For NADP-dependent enzyme and protein
	Activated Affinity (Ligand Immobilization)	TSKgel Tresyl-5PW	TOYOPEARL AF-Tresyl-650M, AF-Formyl-650M TOYOPEARL AF-Epoxy-650M, AF-Amino-650M TOYOPEARL AF-Carboxy-650M	•Tresyl resin requires no additive for reaction
Reversed-Phase Chromatography (RPC)	Protein, Peptide	TSKgel Protein C4-300 TSKgel ODS-120H, ODS-100V TSKgel Octadecyl-NPR, Octadecyl-2PW, -4PW TSKgel Phenyl-5PW RP	N/A	•Separation of protein and peptide with molecular weight below 50,000. •TSKgel Phenyl-5PW RP shows better resolution and recovery for protein with molecular weight over 30,000.
Normal-Phase Chromatography (NPC, HILIC)	Peptide Oligosaccharide, Polysaccharide	TSKgel Amide-80 TSKgel NH2-100	N/A	•Glycan analysis depleted from antibody and glycoprotein





# Attention

- Most of Tosoh products are not inspected for safety and toxicity. Tosoh products should not be assured as harmless nor non-toxic without any specific warning and/or caution.
- When products or those solution purified by Tosoh products are applied to finished or intermediate products, please fully confirm its safety by yourselves.
- Data described on this presentation are Tosoh data or reference data, but not assured ones.
- Reconfirmation or taking data by yourselves should be highly appreciated according to your environment, conditions and judgement criteria.





# Thank you for your attention.

**Web Site;**

<https://www.separations.us.tosohbioscience.com/>

<https://www.separations.eu.tosohbioscience.com/>





# Flow-Through Purification of mAb Fraction on TOYOPEARL NH<sub>2</sub>-750F

Purification Process	Conditions	Monomer Yield (%)	Aggregate (%)*	HCP Log Reduction	DNA (pg/mg)	Protein A (ppm)
TOYOPEARL NH <sub>2</sub> -750F (Flow through)	Start **	N/A	3.1	1.7	0.3	<0.05
	A	95	2.6	2.3	not detected***	<0.05
	B	88	1.0	2.5	not detected***	<0.05

\* Aggregate determined according to Bond et al.

\*\* Starting material; TOYOPEARL AF-rProtein A HC-650F capture eluate with 0.1 mol/L glycine

Sample materials (Protein A eluate) were adjusted with Tris to keep pH 6.8 – 7.0, after 30 min incubation for viral inactivation.

\*\*\* Less than the lower limit of quantification using resDNASEQ Quantitative CHO DNA Kit (Life Tech)

- Improvement of removal of impurities from Protein A eluent by flow-through
- 2.5 Log reduction of HCP
- Host DNA not detected, Protein A leacheable less than 0.05 ppm


Reference; A. Grabski et al., PREP 2016, Poster, Presentation data was slightly modified.

- Flow-through conditions would be more effective at NaCl over 0.15 mol/L in buffer at neutral pH on TOYOPEARL NH<sub>2</sub>-750F

Reference; Tosoh Bioscience, Chromatography Application Note




# Purification of mAb Digests on Ca<sup>++</sup>Pure-HA




Resin: CaPure-HA  
Column size: 6.6 mm ID × 15 cm (5.1 mL)  
Mobile phase: A: 5 mmol/L sodium phosphate, 0.5 mmol/L CaCl<sub>2</sub>, pH 6.5  
B: mobile phase A + 500 mmol/L sodium phosphate  
Gradient: 12% B (5CV), 22.5% B (5CV), 40% B (5CV)  
Flow rate: 2.57 mL/min (450 cm/hr)  
Detection: UV @ 280 nm, conductivity  
Sample: papain digested mAb, 5.0 g/L  
Sample load: 2.0 g/L-resin

# Flow-Through (FT) HIC Purification of mAb (resin under development)

## • Condition-1 : 50 mM Phosphate buffer (5 mS/cm)



**Fig. 2. Chromatogram of Flow through purification**



**Fig. 3. SEC chromatograms of Load and Flow through pool**

Column : TSKgel® G3000SWxL  
(7.8 mm I.D. x 30 cm)  
Mobile phase : 100 mM Phosphate buffer + 200 mM Sodium sulfate pH 6.8  
Flow rate : 1.0 mL/min

**Table 1. Purity and recovery rate by SEC analysis**

	Monomer [%]	Dimer [%]	HMW [%]	Recovery [%]
Load	87.7	6.1	6.2	-
Flow through pool	99.4	0.6	0.0	96

Resin : Prototype HIC resin  
FT buffer : 50 mM Phosphate buffer, 5 mS/cm  
Strip buffer : Acetate buffer, 0.2 mS/cm  
Flow rate : 0.8 mL/min (Residence time = 2.0 min)  
Sample : 5 g/L Humanized mAb @ FT buffer  
Load Vol. : 20 mL