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
# Duolink® PLA Control Kit – PPI (#DUO92202)

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- A. 兩片 8-well chamber slide with EGF-treated pre-fixed SKOV3 cells
- B. Mouse anti-EGFR & rabbit anti-HER2 antibodies

### Duolink® PLA Control Kit – PPI slide layout and sample placement.

All wells contain EGF-stimulated, pre-fixed SK-OV3 cells. After permeabilization and blocking, wells 1 and 5 were incubated with both rabbit anti-EGFR and mouse anti-HER2 primary antibodies, wells 2 and 6 were incubated with anti-EGFR antibody alone, wells 3 and 7 were incubated with anti-HER2 antibody alone, and wells 4 and 8 were incubated with antibody diluent only.

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
Duolink <i>in situ</i> PLA Probes ( 30rxn/100rxn )	
Anti-Mouse PLUS	DUO92001
Anti-Rabbit MINUS	DUO92005
Anti-Goat MINUS	DUO92006
Anti-Rabbit PLUS	DUO92002
Anti-MOUSE MINUS	DUO92004
Anti-Human PLUS	DOU92020
Anti-Human MINUS	DUO92021

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


#### Preparation

Step 1. Prepare microscope slides or microplate


Step 2. Fix cells or tissue, wash, permeabilize, and block


Step 3. Generate PLA probes using Multicolor Probemaker kits


#### Duolink® PLA


Step 4. Add Multicolor PLA probes

Step 5. Wash and add the ligation solution

Step 6. Wash and add the amplification solution


Step 7. Wash and add the detection solution


#### Analysis

Step 8. Capture images of each protein event separately

Step 9. Obtain objective quantification of each protein event separately

Step 10. Combine data analysis

## ➤ Duolink® PLA Multicolor Probemaker Kit (標 20ug x 2 抗體)

Product #	Fluor Color	Ex. (nm)	Em. (nm)	Filter
DUO96010	Red	593	622	Texas Red
DUO96020	Green	490	520	FITC
DUO96030	Orange	542	562	Cy™3
DUO96040	Far Red	646	664	Cy™5

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Product #	Product Name	PKG
DUO96000	Duolink® PLA Multicolor Reagent Pack	30 RXN
	(內含四色螢光, 最多可同時辨認偵測 for 四色螢光之 probe)	100 RXN
DUO82049	Duolink® In Situ Wash Buffer, Fluorescence	4L
		20L
DUO82040	Duolink® In Situ Mounting Medium with DAPI	5ML

## PLA® Validated Antibody Targets

		Target			
14-3-3 Sigma	Chk2	GbetaL	Menin	Rad17	SMN1
53BP1	c-Jun	GNL1	MLH1	Rad17, Phospho (S645)	SMRT
53BP1, Phospho (S25)	CKI alpha	GRLF1/p190RhoGAP	MLL1	Rad18	Sp1
AATF/Che-1	CKII alpha	GTF2B/TFIIB	MOF/MYST1	Rad21	SRC1
AKT1	CKII beta	GTF3C1/TFIIIC220	MOV10	Rad50	SRC2
AKT2	Claspin	H2AX	Mre11	Rad9	SRC3
ARID4B	c-myc	H2AX	MSH2	RanBP3	STAT3
ASH2	Coilin	HBO	MSH6	RbBP4	STUB1/CHIP
ATAD5/FRAG1	CoREST	HCF1	MLH1	RbBP7	TAF5
ATF2	CPSF160	HDAC1	MLL1	RCC1	TCP1
ATM	CREB	HdmX/MDM4	MOF/MYST1	RCD8	TDP1
ATMIN	c-Rel	HdmX/MDM4, Phospho (S403)	MOV10	RelA	TIF1 Alpha/TRIM24
ATR	CRM1	HIF1-alpha	Mre11	RelA, Phospho (S468)	Timeless
ATRIP	CRSP1/TRAP220	HMG20B/BRAF35	MSH2	RelB	TIPIN
ATRX	CSN1	HP1-beta	MSH6	RENT1	TLK1
BAF53A	CTBP1	HSP60	mTOR	RFC1	TLK2
BAP1	CTF18	IKK-alpha	NBS1	RENT1	TNIK
BARD1	CtIP	IKK-beta	NCoR	RFC1	TopBP1
BCAR3	Cul1	IQGAP1	NEK9	RFC2	TRF2
BCR	Cul4a	JARID1A/RBP2	NF-kappaB1	RFC3	TRRAP
BIRC6/Apollon	Cyclin T1	KAP-1	NF-kappaB2	RFC4	TRX2
BLAP75	DAXX	KAP-1, Phospho (S824)	NFRKB	RFC5	TSG101
BLM	DBC1/p30 DBC	Ki-67	NIF1	Rictor	TTI1
BMI1	DBF4	KPNA3	NONO	RIF1	USP28
BRCA1	DC8	Ku70	NOP56	RNA Polymerase II	USP7
BRCA1, Phospho (S1524)	Dcp2	Ku80	Notch1	RNA Polymerase II, Phospho (S2)	VCP
BRD2	DDB1	Lasu1/Urb1	NPM1	ROCK2	Vimentin
BRD4	DDX56	LATS2	ORC1	RPA32	WAPL
BRF1	DMAP1	MAD1	ORC2	RPA32, Phospho (S33)	WDR5
BRG1/SMARCA4	DNA-PKcs	MAD2	p16INK4a	RPA32, Phospho (S4/S8)	XRN1
cAbl	DNMT1	MCM2	p300	RPA70	YL1
Caf1/p150	E2F1	MCM2, Phospho (S108)	p400	RPL26	YY1
CBP	E2F4	MCM2, Phospho (S27)	p53	RPS6, Phospho (S235/S236)	ZNF8
CBX3	EGFR	MCM2, Phospho (S40/S41)	p73	RuvBL2	—
CBX5	EPC1	MCM2, Phospho (S41)	PARP1	SA1	—
CCDC131	ERK1	MCM2, Phospho (S53)	PCNA	SCC-112	—
CDC25a	ERK2	MCM3	PHF8	SDS3	—
CDC25c	ESCO2	MCM4	PIM1	SETDB1	—
CDC6	Estrogen Receptor Alpha	MCM3	PKM2	Sin1	—
CDC7	Estrogen Receptor Alpha, Phospho (S305)	MCM4	PLK1	Sin3A	—
CDK1	EWS	MCM5	PML	Sirt1	—
CDK2	FANCD2	MCM6	PMS2	SMC1	—
CDK7	FANCE	MCM7	PoIA1	SMC1, Phospho (S957)	—
CDK9	FOXO1a	MCPH1	PPP1CA	SMC1, Phospho (S966)	—
CHD4/Mi2 beta	FOXO3a	MDC1	PPP1CC	SMC2	—
Chk1	FUS	MEK1	PPP5C	SMC3	—
Chk1	G3BP1	MEK1, Phospho (T286)	PRKCBP1	SMC4	—
Chk1, Phospho (S317)	GAPDH	MEK1	PSF/SFPQ	SMG1	—

Enabling Science to  
Improve the Quality of Life

Order/Customer Service: [sigma-aldrich.com/order](http://sigma-aldrich.com/order)  
 Technical Service: [sigma-aldrich.com/techservice](http://sigma-aldrich.com/techservice)  
 Development/Custom Manufacturing Inquiries **SAFC**® [safcglobal@sial.com](mailto:safcglobal@sial.com)  
 Safety-related Information: [sigma-aldrich.com/safetycenter](http://sigma-aldrich.com/safetycenter)

World Headquarters  
 3050 Spruce St.  
 St. Louis, MO 63103  
 (314) 771-5765  
[sigma-aldrich.com](http://sigma-aldrich.com)

©2014 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA, SAFC and SIGMA-ALDRICH are trademarks of Sigma-Aldrich Co. LLC, registered in the US and other countries. Where bio begins is a trademark of Sigma-Aldrich Co. LLC. Duolink and PLA are registered trademarks of Olink AB. Sigma brand products are sold by affiliated Sigma-Aldrich distributors. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see product information on the Sigma-Aldrich website at [www.sigmaaldrich.com](http://www.sigmaaldrich.com) and/or on the reverse side of the invoice or packing slip.