

3D Cell Culture 天然的水凝膠



購買水凝膠滿1萬再送
血清瓶造型隨身碟1支

水凝膠常見各種實驗應用


Tumor Spheroids



Organoids Culture



Invasion assay



Tube formation assay



ES & iPSC Culture



ACD kit for Tumor Spheroid


- 天然動物萃取成分，專為Tumor Spheroid 形成所設計配方
- 獨特的4度孵育才成膠，免除使用Matrigel需要預冷所有材料的不便性
- 搭配套組內的D buffer, 讓您輕鬆就能回收Tumor Spheroid 超簡單




24 assay **\$4500**
XR-LS01A01-1KIT

480 assay **\$15000**
XR-LS01A02-1KIT

Gelation Protocol:




Dissolving Protocol:



適用的細胞株

Cell lines origin	Name
Tongue (human)	SCC-15
Laryngeal (human)	Hep-2
Muscle (human)	RD
Breast (human)	MCF7, ZR-75-1, MDA-MB-231, MDA-MB-361
Lung (human)	H292, A549
Liver (human)	Huh7, Hep G2
Kidney (human)	CRL-1932
Prostate (human)	DU145
Stomach (human)	AGS, SCM-1
Pancreas (human)	PANC-1, MIA PaCa-2, ASPC-1
Colon (myofibroblasts, human)	CT5.3
Colorectal (human)	LoVo, SW620, HT-29, HCT 116
Ovary (human)	SKOV3, OVCAR-3
Cervix (human)	NHIK 3025

Huh7



A549

