

高效核酸萃取套組

 日期 **3/6~4/30**


- ✓ 大大提升核酸純度、濃度與產量
- ✓ 增量30%檢體量，可處理高達 1000 μ l
- ✓ 檢體洗脫體積 $\leq 30 \mu$ l，樣品濃度提高 **3** 倍
- ✓ 簡化操作步驟，減少離心次數
- ✓ 應用廣泛：萃取出核酸適用於後續的 PCR、qPCR、次世代定序



品名	品號	次數
FavorPrep™ Tissue Genomic DNA Extraction HE Mini Kit	FATG1033	
FavorPrep™ Blood Genomic DNA Extraction HE Mini Kit	FABG1033	
FavorPrep™ Stool DNA Extraction HE Mini Kit	FAST1033	
FavorPrep™ Soil DNA Extraction HE Mini Kit	FASO1033	50preps
FavorPrep™ Urine DNA Extraction HE Mini Kit	FAUD1033	
FavorPrep™ Bacterial Genomic DNA Extraction HE Mini Kit	FABD1033	
FavorPrep™ Bacterial RNA Extraction HE Mini Kit	FABR3033	
FavorPrep™ Circulating Nucleic Acid Extraction HE Mini Kit	FACN1033	



 1.5ml 反應管
 #SAR-72.706.201-1CS


- ✓ 單管獨立包裝適合NGS送樣
- ✓ 最高等級純淨度認證
- ✓ 無核酸與核酸酶
- ✓ 無致熱源
- ✓ ISO11135滅菌驗證



- Human DNA <0.5 fg/ μ l
- Bacterial DNA <0.2 fg/ μ l
- DNase <5x10⁻⁷U / μ l
- RNase <5x10⁻¹¹ Kunitz units/ μ l
- ATP <1x10⁻¹² mmol/ μ l
- Pyrogens <0.002 EU/ml
- Sterility validated (ISO 11135)

 特價 **1200**/CS(120支) 數量有限售完為止

<p>特價 3000 500/CS</p>  <p>滅菌</p> <p>13ml 養菌管 SAR-62.515.006-1CS</p>	<p>特價 1890 300/CS</p>  <p>滅菌 低內毒素</p> <p>50ml 離心管 SAR-62.547.254-1CS</p>	<p>特價 2100 300/CS</p>  <p>滅菌 低內毒素</p> <p>50ml 可站立離心管 SAR-62.559.001-1CS</p>	<p>特價 2300 500/CS</p>  <p>滅菌 低內毒素</p> <p>15ml 離心管 SAR-62.554.502-1CS</p>	<p>特價 2000 100支</p>  <p>滅菌</p> <p>30ml 螺旋管 SAR-62.543.001-BAG</p>
<p>特價 1000</p>  <p>20孔試管架 (適用 30ml 螺旋管) SAR-93.841.100-1EA</p>	<p>特價 9000 500/CS</p>  <p>唾液收集管 附棉花棒 SAR-51.1534-1CS</p>	<p>特價 5500 1000/CS</p>  <p>5ml 螺旋管 SAR-60.558-1CS</p>	<p>特價 1000 500/BAG</p>  <p>螺旋管蓋 SAR-65.716-1BAG</p>	<p>特價 1500 500/BAG</p>  <p>滅菌</p> <p>螺旋管蓋 SAR-65.716.723-1BAG</p>
<p>特價 10000 5000/CS</p>  <p>1.5ml 螺旋管 SAR-72.703-1CS</p>	<p>特價 5000 2500/CS</p>  <p>滅菌</p> <p>2ml 螺旋管 SAR-72.609.711-1CS</p>	<p>特價 10000 5000/CS</p>  <p>1.5ml 螺旋管 SAR-72.692-1CS</p>	<p>特價 5800 1000/CS</p>  <p>滅菌</p> <p>2ml 螺旋管 SAR-72.693.105-1CS</p>	<p>特價 6000 1000/CS</p>  <p>滅菌</p> <p>2ml 螺旋管 SAR-72.694.006-1CS</p>
<p>特價 1200 120/CS</p>  <p>滅菌</p> <p>1.5 ml 反應管 SAR-72.706.201-1CS</p>	<p>特價 5000 1000/CS</p>  <p>滅菌</p> <p>0.5ml 螺旋管 SAR-72.730.005-1CS</p>	<p>特價 5500 1000/CS</p>  <p>滅菌</p> <p>0.5ml 螺旋管 SAR-72.730.006-1CS</p>	<p>特價 5500 500/CS</p>  <p>滅菌</p> <p>0.5ml 螺旋管 SAR-72.730.217-1CS</p>	<p>特價 4000 24/CS</p> <p>v底設計 2.2ml</p>  <p>96 well 深孔盤 SAR-82.1972-1CS</p>
<p>特價 2500 1000/CS</p>  <p>滅菌</p> <p>10ul 接種環 SAR-86.1562.010-1CS</p>	<p>特價 2500 1000/CS</p>  <p>滅菌</p> <p>1ul 接種環 SAR-86.1567.010-1CS</p>	<p>特價 2500 1000/CS</p>  <p>滅菌</p> <p>接種針 SAR-86.1568.010-1CS</p>	<p>5個 特價 6500</p> <ul style="list-style-type: none"> 可高溫高壓滅菌 可以重複使用  <p>Silicon insert, 8 well SAR-94.6032.039-1PK</p>	<p>5個 特價 6000</p> <ul style="list-style-type: none"> 可高溫高壓滅菌 可以重複使用  <p>Silicon insert, 4 well SAR-94.6034.067-1PK</p>
<p>特價 800 100/Box</p>  <p>體積: 4.2ml, 材質: PMMA optical sides: 2</p> <p>比色管(雙面透光) SAR-67.738-1BOX</p>	<p>特價 500 100/Box</p>  <p>體積: 3 ml 材質: PS optical sides: 2</p> <p>比色管(雙面透光) SAR-67.742-1BOX</p>	<p>特價 700 100/Box</p>  <p>體積: 4.2ml 材質: PS optical sides: 4 用途: 螢光分析</p> <p>比色管(4面透光) SAR-67.754-1 BOX</p>	<p>特價 800 100/Box</p>  <p>體積: 4.2ml 材質: PMMA optical sides: 4 用途: 螢光分析</p> <p>比色管(4面透光) SAR-67.755-1BOX</p>	<p>特價 1200 100/Box</p>  <p>體積: 2.7ml 材質: 特殊塑料 optical sides: 2 用途: UV應用</p> <p>UV 比色管 2.7ml SAR-67.758-1BOX</p>

DM11503004



友和貿易股份有限公司
UNI-ONWARD Corp.
客服專線: (02)2600-0611

◆台北 TEL:(02)2600-6699 FAX:(02)2600-0799
◆新竹 TEL:(03)658-0508 FAX:(03)658-0509
◆台中 TEL:(04)2355-2008 FAX:(04)2355-2006
◆嘉義 TEL:(05)277-2262 FAX:(05)277-2263

◆台南 TEL:(06)249-5488 FAX:(06)249-5499
◆高雄 TEL:(07)374-8686 FAX:(07)374-8585
◆花蓮 TEL:(03)836-0065 FAX:(03)836-0075
www.uni-onward.com.tw

