


GoPAGE Precast Gel 產品選擇指南



| TGN 系列 (用Tris-Glycin buffer) | 品號 | 尺寸 | 適用電泳槽 | 包裝 |
|-------------------------------------|--------------------|--------------|--------------------------------------|------|
| GoPAGE TGN Precast gel (10%) | XR-LSB-GL4210-1SET | 10X8.3 cm | Bio-Rad | 10/盒 |
| GoPAGE TGN Precast gel (4-15%) | XR-LSB-GL4510-1SET | | | 10/盒 |
| GoPAGE TGN Precast gel (10%) | XR-LSB-GL5210-1SET | 10X10 cm | Novex system, Hoefler SE260 | 10/盒 |
| GoPAGE TGN Precast gel (4-15%) | XR-LSB-GL5510-1SET | | | 10/盒 |
| Bis-Tris 系列 (用MOPS buffer) | 品號 | 尺寸 | 適用電泳槽 | 包裝 |
| GoPAGE Bis-Tris Precast gel (8%) | XR-LSB-GL2110-1SET | 10X8.3 cm | Bio-Rad | 10/盒 |
| GoPAGE Bis-Tris Precast gel (12%) | XR-LSB-GL2310-1SET | | | 10/盒 |
| GoPAGE Bis-Tris Precast gel (4-12%) | XR-LSB-GL2510-1SET | | | 10/盒 |
| GoPAGE Bis-Tris Precast gel (8%) | XR-LSB-GL3110-1SET | 10X10 cm | Novex System, Hoefler SE260 | 10/盒 |
| GoPAGE Bis-Tris Precast gel (12%) | XR-LSB-GL3310-1SET | | | 10/盒 |
| GoPAGE Bis-Tris Precast gel (4-12%) | XR-LSB-GL3510-1SET | | | 10/盒 |